

THE SPANNER™ TEMPORARY PROSTATIC STENT
ORDERING INFORMATION

TO PLACE AN ORDER:

TELEPHONE: Call Customer Service at 888-528-9073

FAX: 218-338-6710

EMAIL: customerservice@abbeymoormedical.com

MAIL: AbbeyMoor Medical, Inc.
Attn: Customer Service
501 East Soo Street
Parkers Prairie, MN 56361
USA

CONFIRMATION: Confirmation of order will be faxed or emailed back to you upon receipt of the order.

SHIPPING AND VAT: Paid by customer where applicable. If you would like AbbeyMoor Medical to use your shipping agent please supply the appropriate information.

CUSTOMER INFORMATION:

Institution: _____ Contact Name: _____
Address: _____
City: _____ Country: _____ Postal Code: _____
Telephone: _____ Fax: _____ Email: _____
Customer VAT number: _____

CREDIT CARD INFORMATION:

Credit Card – (please check): Visa MasterCard American Express

Credit Card Number: _____ Expiration Date: _____

Name of Card Holder: _____

ORDERING INFORMATION:

DESCRIPTION	SPANNER SIZE	QUANTITY
Spanner	5	
Spanner	6	
Spanner	7	
Spanner	8	
Spanner	9	
Note: Each Spanner includes a Surveyor		

RECOMMENDED INITIAL ORDER:

SPANNER SIZE	QUANTITY
5	1
6	1
7	1
8	1

TRUS Measurement (Prostatic Urethra Length— Distance from Bladder Neck to Prostatic Apex)	Spanner Size
2.5 - 3.0	5
3.1 - 3.6	6
3.7 - 4.6	7
4.7 - 5.6	8
5.7 - 6.6	9

